


Performance Characteristics

Uniform coverage	95 % amplicons detected $\geq 0.2x$ mean coverage
High target specificity	96 % on target read counts



Graph presenting the high uniform coverage for all PED MASTR Plus amplicons. To allow comparison between samples, the read pair counts were normalized.

Workflow

SNV and CNV* variant calling: advised maximum number of samples per run

Illumina MiSeq®	
Reagent Kit	Nr samples
Kit v2	5
Kit v3	10



Order info

Cat. No.	Product Name	Samples
MR-0260.024	PED MASTR Plus	24

MID (Molecular Identifiers) kits are necessary to complete the workflow.

© Agilent Technologies, Inc. 2017
 Printed in Belgium, December 31, 2017
 5991-8389EN
 PR7000-1414

For Research Use Only. Not for use in diagnostic procedures.